


Attends® Flex



Attends - The system with 10 absorption levels for bladder weakness

Attends Flex are a range of all-in-one belted pads designed to manage heavy urinary and/or faecal incontinence.



- | | | |
|---|--|--|
| <p>1 White breathable back sheet to aid containment and for improved comfort, discretion and skin health.</p> | <p>2 Wetness indicator aids good nursing practice and acts as a guide to when the pad requires changing.</p> | <p>3 Elasticated waist front and back for improved fit and helps to reduce the risk of leakage.</p> |
| <p>4 Flexible fixation system allows the belt to be fitted snugly around the waist before fitting the pad for ease of fitting.</p> | <p>5 Curved leg elastics for improved fit and to help reduce the risk of leakage.</p> | <p>6 Full barrier leg cuff reduces the risk of leakage.</p> |
| <p>7 Anatomically shaped core for a close comfortable fit.</p> | <p>8 Superabsorbent core with acquisition layer for improved leakage protection, skin dryness and odour protection.</p> | <p>9 The size and absorbency level are printed on the back sheet for easy identification of products in care homes/wards.</p> |

For more information visit www.attends.co.uk

Absorption / Size

Absorption / Size		Packcount Casecount	EAN 7332152 Pack / Case	Hip / Waist Size (cm)	Backsheet	Absorption Core Material	Topsheet Material	Elastics	Leakage Barriers	Resealable fixation system	Natural Odour Protection	Breathable	Quick-Dry Technology
8	Attends Flex Small 	28 4	206713 206720	60 - 85	Laminate of nonwoven and breathable film	Dual Core system with fluff and superabsorbents with Quick-Dry acquisition layer	Non-woven	✓	✓	✓	✓	✓	✓
	Attends Flex Medium 	28 4	206775 206782	75 - 105	Laminate of nonwoven and breathable film	Dual Core system with fluff and superabsorbents with Quick-Dry acquisition layer	Non-woven	✓	✓	✓	✓	✓	✓
	Attends Flex Large 	28 4	206836 206843	85 - 125	Laminate of nonwoven and breathable film	Dual Core system with fluff and superabsorbents with Quick-Dry acquisition layer	Non-woven	✓	✓	✓	✓	✓	✓
	Attends Flex Extra Large 	26 4	206898 206904	100 - 155	Laminate of nonwoven and breathable film	Dual Core system with fluff and superabsorbents with Quick-Dry acquisition layer	Non-woven	✓	✓	✓	✓	✓	✓
9	Attends Flex Small 	26 4	206737 206744	60 - 85	Laminate of nonwoven and breathable film	Dual Core system with fluff and superabsorbents with Quick-Dry acquisition layer	Non-woven	✓	✓	✓	✓	✓	✓
	Attends Flex Medium 	26 4	206799 206805	75 - 105	Laminate of nonwoven and breathable film	Dual Core system with fluff and superabsorbents with Quick-Dry acquisition layer	Non-woven	✓	✓	✓	✓	✓	✓
	Attends Flex Large 	26 4	206850 206867	85 - 125	Laminate of nonwoven and breathable film	Dual Core system with fluff and superabsorbents with Quick-Dry acquisition layer	Non-woven	✓	✓	✓	✓	✓	✓
	Attends Flex Extra Large 	22 4	206911 206928	100 - 155	Laminate of nonwoven and breathable film	Dual Core system with fluff and superabsorbents with Quick-Dry acquisition layer	Non-woven	✓	✓	✓	✓	✓	✓
10	Attends Flex Small 	22 4	206751 206768	60 - 85	Laminate of nonwoven and breathable film	Dual Core system with fluff and superabsorbents with Quick-Dry acquisition layer	Non-woven	✓	✓	✓	✓	✓	✓
	Attends Flex Medium 	22 4	206812 206829	75 - 105	Laminate of nonwoven and breathable film	Dual Core system with fluff and superabsorbents with Quick-Dry acquisition layer	Non-woven	✓	✓	✓	✓	✓	✓
	Attends Flex Large 	22 2	206874 208007	85 - 125	Laminate of nonwoven and breathable film	Dual Core system with fluff and superabsorbents with Quick-Dry acquisition layer	Non-woven	✓	✓	✓	✓	✓	✓
	Attends Flex Extra Large 	20 4	206935 206942	100 - 155	Laminate of nonwoven and breathable film	Dual Core system with fluff and superabsorbents with Quick-Dry acquisition layer	Non-woven	✓	✓	✓	✓	✓	✓



Quick-Dry

Quick Dry is the acquisition layer placed directly below the top sheet where the urine enters the pad. It quickly absorbs urine away from the skin and into the pad and prevents the urine returning to the surface. It can handle several gushes of urine without affecting how the pad works which results in exceptional skin dryness and leakage protection for the user.



Active Zone

The Active Zone is an area centred in the middle of the pad, which provides high absorption where it is really needed. This provides exceptional skin dryness and improved leakage protection for the user.



Natural Odour Protection

Natural Odour Protection reduces the risk of unpleasant smells.



pH Value

The core technology with the curly fibre acquisition layer helps to guarantee a low pH which is similar to the skin. This helps to protect the skin and ensures improved skin health.



Breathable

The breathable material enables the skin to breathe, which contributes to improved skin health.



Attends Flex are tested and approved by proDERM, Institute for Applied Dermatological Research.