


Attends® Slip Active



Attends - The system with 10 absorption levels for bladder weakness

Attends Slip Active are a range of discreet all-in-one pads designed to manage heavy urinary and/or faecal incontinence.



- | | | |
|--|--|---|
| 1 White waterproof back sheet to aid containment. | 2 Wetness indicator aids good nursing practice and acts as a guide to when the pad requires changing. | 3 Breathable textile like wings provide comfort and support around the waist as well as improved skin health. |
| 4 Re-sealable two step fixation system for ease of fitting. Use the blue positioning tape to preapply the product and adjust if necessary. Then fasten the wide tapes to secure the pad in place. | 5 Curved leg elastics for improved fit and to help reduce the risk of leakage. | 6 Full barrier leg cuff reduces the risk of leakage. |
| 7 Anatomically shaped core for a close comfortable fit. | 8 Super absorbent core with acquisition layer for improved leakage protection, skin dryness and odour protection. | 9 The size and absorbency level are printed on the back sheet for easy identification of products in care homes/wards. |

For more information visit www.attends.co.uk

Absorption / Size

Absorption / Size		Packcount Casecount	EAN 7332152 Pack / Case	Hip / Waist Size (cm)	Backsheet	Absorption Core Material	Topsheet Material	Breathable	Elastics	Leakage Barriers	Resealable fixation system	Natural Odour Protection	Quick-Dry Technology
8	Attends Slip Active Medium 	28 4	200728 200735	75 - 110	Polyethylene film	Dual Core system with fluff and superabsorbents with Quick-Dry acquisition layer	Non-woven	Yes - Wings	✓	✓	✓	✓	✓
	Attends Slip Active Large 	28 4	201329 201336	100 - 145	Polyethylene film	Dual Core system with fluff and superabsorbents with Quick-Dry acquisition layer	Non-woven	Yes - Wings	✓	✓	✓	✓	✓
9	Attends Slip Active Medium 	28 4	201503 201510	75 - 110	Polyethylene film	Dual Core system with fluff and superabsorbents with Quick-Dry acquisition layer	Non-woven	Yes - Wings	✓	✓	✓	✓	✓
	Attends Slip Active Large 	28 4	201381 201398	100 - 145	Polyethylene film	Dual Core system with fluff and superabsorbents with Quick-Dry acquisition layer	Non-woven	Yes - Wings	✓	✓	✓	✓	✓
10	Attends Slip Active Medium 	28 4	201404 201411	75 - 110	Polyethylene film	Dual Core system with fluff and superabsorbents with Quick-Dry acquisition layer	Non-woven	Yes - Wings	✓	✓	✓	✓	✓
	Attends Slip Active Large 	28 2	201428 201435	100 - 145	Polyethylene film	Dual Core system with fluff and superabsorbents with Quick-Dry acquisition layer	Non-woven	Yes - Wings	✓	✓	✓	✓	✓



Quick-Dry

Quick-Dry is the acquisition layer placed directly below the top sheet where the urine enters the pad. It quickly absorbs urine away from the skin and into the pad and prevents the urine returning to the surface. It can handle several gushes of urine without affecting how the pad works which results in exceptional skin dryness and leakage protection for the user.



pH Value

The core technology with the curly fibre acquisition layer helps to guarantee a low pH which is similar to the skin. This helps to protect the skin and ensures improved skin health.



Natural Odour Protection

Natural Odour Protection reduces the risk of unpleasant smells.



Active Zone

The Active Zone is an area centered in the middle of the pad, which provides fast absorption where it is really needed. This provides exceptional skin dryness and improved leakage protection for the user.



Attends Slips are tested and approved by proDERM, Institute for Applied Dermatological Research.