

0 1 2 3 4 5 6 7 8 9 10

The Attends numbering system - Absorbency levels

Question:

When should I use a washable product?

What is special about Attends Washable products?

Answer:

Attends Washables are a range of products designed to manage light to moderate urinary incontinence.

The products come in an range of sizes and absorbencies. The range also includes a bed pad range.




Attends washable products comprise of male, female and unisex absorbent brief, a children's range for boys and girls which are available in different sizes to meet client needs. The range also includes two absorbent bed pads.

The briefs look like normal underwear and have a built in absorbent pad. This allows the pad to stay in place and prevents twisting of the pad especially during the washing cycle.



The product has a waterproof vinyl backing which acts as a barrier and prevents wetness leaking on to clients clothing or bedding. A soaker layer draws urine away from the top sheet by absorbing and wicking to give improved dryness and reduce the risk of leakage. In addition MicroSafe (apart from the 2 litre bed pad), an antimicrobial fibre effectively, defends against odour and infection providing the client with greater confidence and independence.

The products are guaranteed for 300 washes and can be washed at temperatures up to 95°C. The DSS recommends that any wash over 76°C will kill bacteria and control infection.

Attends Washable Brief Adults For Women: the measurements apply to hip size. For Men: the measurements apply to waist size

	Absorbency	Small	Medium	Large	Extra Large	Extra Extra Large
Size		81-86cm	91-96cm	101-106cm	111-117cm	122-132cm
Male 	200mls	20 18 55	20 18 79	20 18 93	20 19 16	20 19 30
Female 	270mls	20 19 54	20 19 78	20 19 92	20 20 12	20 20 36
Unisex 	500 mls	20 20 50	20 20 74	20 20 98	20 21 11	20 21 35
Pack Count		1	1	1	1	1

Attends Washable Brief Children For Children: the measurements apply to waist size

	Absorbency	Age 3-4yrs	Age 5-6yrs	Age 7-8yrs	Age 9-10yrs	Age 11-12yrs
Size		46-56cm	51-61cm	56-66cm	61-71cm	71-81cm
Girls 	150 mls	20 21 59	20 21 73	20 21 97	20 22 10	20 22 34
Boys 	150mls	20 22 58	20 22 72	20 22 96	20 23 19	20 23 33
Pack Count		1	1	1	1	1

Attends Washable Bed Pad

Size	2 Litre	3 Litre
7332152	20 23 57	20 23 71
Pack Count	1	1

For more information visit www.attends.co.uk