

Attends

For Men

0	1	2	3	4	5	6	7	8	9	10
---	---	---	---	---	---	---	---	---	---	----

The Attends numbering system – Absorption degree

Question:

What happens?

Answer:



Loss of a few drops of urine after the main stream, when the bladder appears to be empty. Many men relate to the term 'after dribble' or 'post micturition dribbling'

Attends for men is a new anatomically shaped drip shield and is available in two absorption degrees. A "tighter" core, named airlaid, provides a product that is slim and yet highly absorbent.

Odour Protection technology minimizes the risk of unpleasant odours developing.

Attends for men has a cottonlike material on the backsheet, giving the feeling of a natural addition to underwear. The adhesive strips keep the product in place.



	
Attends for Men	Attends for Men
I	II

Art. nr.-case		899561	201794
Length	mm	190	190
Thickness	mm	1.5	2.6
Weight	g	6	10
ISO Absorption (acc. to ISO 11948)	g	110	170
Topsheet	Mtrl.	Non-Woven	Non-Woven
Core		Airlaid-core incl. SAP	Airlaid-core incl. SAP
Backsheet	Mtrl.	Polyethylene with nonwoven	Polyethylene with nonwoven
Odour Protection		yes	yes
Adhesive Strips		yes	yes
Packcount	Pcs	25	16
Cascount	Packs	4	4
Nr of cases/pallet		80	80