

Attends® For Men





Attends - The system with 10 absorption levels for bladder weakness

Attends For Men are anatomically shaped shields designed for men with urinary incontinence.

Levels 1 and 2 are designed for men with post micturition dribble.

Levels 3 and 4 for men with light to moderate urinary incontinence.







- 1 Anatomically shaped for a close comfortable fit for the male anatomy and improved leakage protection.
- 2 Fixation tapes that allow the product to be worn with close fitting underwear and help to reduce the risk of leakage by ensuring the pad stays securely in place.
- 3 A soft breathable textile backsheet reduces the risk of rustle for extra discretion and improved comfort.
- 4 Airlaid core enables the product to be very thin for discretion yet highly absorbent for improved leakage protection, skin dryness and odour protection.

- 1 Anatomically shaped for a close comfortable fit for the male anatomy and improved leakage protection.
- 2 Fixation tapes that allow the product to be worn with close fitting underwear and help to reduce the risk of leakage by ensuring the pad stays securely in place.
- 3 A soft breathable textile backsheet reduces the risk of rustle for extra discretion and improved comfort.
- 4 Super absorbent core for improved leakage protection, skin dryness and odour protection.
- 5 Leakage barriers for improved fit and to reduce the risk of leakage.

For more information visit www.attends.co.uk

Attends For Men

	Packcount Casecount	EAN 7332152 Pack / Case	Backsheet	Absorption Core Material	Topsheet Material	Leakage Barriers	Fixation Tapes	Natural Odour Protection	Breathable
Attends For Men 1 	25 4	410032 899561	Laminate of nonwoven and breathable film	Airlaid core including superabsorbents	Non-woven	X	✓	✓	✓
Attends For Men 2 	16 4	201787 201794	Laminate of nonwoven and breathable film	Airlaid core including superabsorbents	Non-woven	X	✓	✓	✓
Attends For Men 3 	14 12	205945 205952	Laminate of nonwoven and breathable film	Fluff and superabsorbents	Non-woven	✓	✓	✓	✓
Attends For Men 4 	14 12	205969 205976	Laminate of nonwoven and breathable film	Fluff and superabsorbents	Non-woven	✓	✓	✓	✓



Natural Odour Protection

Natural Odour Protection reduces the risk of unpleasant smells.



Breathable

The breathable back sheet enables the skin to breathe, which contributes to improved skin health.