


Attends® Slip Regular Plus







Attends - The system with 10 absorption levels for bladder weakness

Attends Slip Regular Plus are a range of all-in-one pads designed to manage heavy/severe urinary and/or faecal incontinence.



- | | | |
|--|--|---|
| 1 White waterproof back sheet to aid containment. | 2 Wetness indicator aids good nursing practice and acts as a guide to when the pad requires changing. | 3 Resealable fixation system with four tapes to secure the pad in place. |
| 4 Elasticated waistband for improved fit. | 5 Curved leg elastics for improved fit and to help reduce the risk of leakage. | 6 Leakage barriers to reduce the risk of leakage. |
| 7 Anatomically shaped core for a close comfortable fit. | 8 Super absorbent core with acquisition layer for improved leakage protection, skin dryness and odour protection. | 9 The size and absorbency level are printed on the back sheet for easy identification of products in care homes/wards. |

For more information visit www.attends.co.uk

	Packcount Casecount	EAN 7332152 Pack / Case	Hip / Waist Size (cm)	Backsheet	Absorption Core Material	Topsheet Material	Elastics	Leakage Barriers	Resealable fixation system	Natural Odour Protection 
Attends Slip Regular Plus Medium 10 	14 4	201442 201459	80 - 115	Polyethylene film	Dual Core system with fluff and super absorbents with acquisition layer	Non-woven	✓	✓	✓	✓
Attends Slip Regular Plus Large 10 	14 4	201466 201473	115 - 150	Polyethylene film	Dual Core system with fluff and super absorbents with acquisition layer	Non-woven	✓	✓	✓	✓
Attends Slip Regular Plus Extra Large 10 	14 4	201480 201497	150 - 175	Polyethylene film	Dual Core system with fluff and super absorbents with acquisition layer	Non-woven	✓	✓	✓	✓



Natural Odour Protection

Natural Odour Protection reduces the risk of unpleasant smells.